

Connect 12 v2.1 – What's New?

Features and Benefits

OS upgrade to Android 12 (Snow Cone)

- Faster boot time
- Run latest apps (educational, professional, social, and personal)
- Supports modern security protocols
- Improvements added from OS6 to OS12, examples:
 - Enhanced accessibility
 - Dark theme
 - Increased battery life
 - Quicker start up
 - Learn more: [Android Releases](#)
- New OS preferred and supported by IT professionals

Prodigi v4.7

- Upgraded to run on Android 12 OS
- Industry leading low vision productivity software
- Delivered on a modern OS with all its new features
- Review for yourself why Connect 12 is the total package:
 - [Connect 12 Webpage](#)
 - [Connect 12 Video](#)

Updated Wi-Fi hardware

- IPv4 and IPv6 compatible
- Access networks using the following:
 - Wi-Fi: 802.11a/b/g/n/ac
 - 1x1 Access Points with SISO standard
- Support for modern Wi-Fi security
- Better connectivity across networks
- Better connectivity for distance cameras
- Preferred and supported by IT professional

Updated Bluetooth hardware

- Support for Bluetooth 5.0
- Better connectivity for peripheral devices
- Preferred and supported by IT professionals

New USB-C power port

- Easy to identify and connect
- Robust, rugged, and sturdy
- Replaceable standard cable

[OPTION] Microsoft Adaptive Kits

- Support persons who need visual and tactile indicators
- Broaden accessibility to other disabilities
- Identify critical buttons and ports
- Use with related peripherals (keyboard, camera, etc.)
- Key components:

- Bump labels
 - Keycap labels and applicator
 - Port labels
 - Opening support
- Free to first 250 FULL v2.1 product purchasers
- Available as an accessory purchase
- Learn more: [Surface Adaptive Kit](#)